# BIO GENE

## HSV PCR

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1. Aldea, C; Alvarez, C.P; Folgueira, L, et al. Detection, quantification and genotyping of Herpes Simplex Virus in cervicovaginal secretions by real-time PCR: a cross sectional survey. *J Clin Microbiol.* 2002 Mar;40(3):1060-2.
2. Chudzik, E; Karabin, K; Dzieciotkowski, T, et al. [Development of a novel real-time PCR method for detection of HSV types 1 and 2 DNA using HybProbe chemistry]. *Med Dosw Mikrobiol.* 2010;62(3):255-62. Polish.
3. CLINICAL AND LABORATORY STANDARDS INSTITUTE (CLSI). Collection, transport, preparation and storage of specimens for molecular methods; approved guideline. CLSI document MM13-A. Pennsylvania, USA: Clinical and Laboratory Standards Institute, 2005.
4. Kessler, H.H; Mühlbauer, G; Rinner, B, et al. Comparison of nested PCR and real time PCR of Herpesvirus infections of central nervous system in HIV patients. *J Clin Microbiol.* 2000 Jul;38(7):2638-42.

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Review: August/2016

# INDEX

Function .....	3
Principle of Action .....	3
Presentation .....	3
Reagents .....	4
Equipments and Operational Inputs .....	4
Transportation and Storage Conditions .....	4
Special Care .....	5
Samples .....	5
Process Description .....	6
A. DNA Extraction .....	6
B . Preparation of the reagents .....	6
C . Dilution of Quantitative Standard .....	7
D . PCR Preparation .....	7
E . Thermocycler settings for Real-Time PCR .....	8
F . Result Validation .....	9
G.Result Interpretation.....	10
Process Limitations .....	11
Product Performance / Quality Control .....	11
Comparison of Methods and Methodology Specificity .....	11
Repeatability .....	11
Reproducibility .....	12
Clinical Sensitivity .....	12
Analytical Sensitivity .....	12
Diagnostic Significance .....	12
Bibliographic References .....	13
Customer Service .....	13
Universal Symbology .....	14

**FUNCTION**

Quantitative detection of *Herpes Simplex Virus* (HSV) types 1 and 2 DNA using real-time Polymerase Chain Reaction (PCR) technology. Designed for *in vitro* diagnostic use only.

**PRINCIPLE OF ACTION**

The **Bio Gene HSV PCR** test is a *in vitro* assay based on the quantitative detection of HSV DNA by real-time PCR.

The Real Time PCR method is used to amplify the pathogen DNA.

A thermocycler for Real Time PCR is used to amplify and detect the fluorescent probe. Software calculates the HSV DNA concentration, expressed in copies/ $\mu\text{L}$ , using the standard curve generated from the quantitative standard, contained in the kit.

Reagent	Presentation
R1	1 x 55 $\mu\text{L}$ *
R2	1 x 525 $\mu\text{L}$ *
R3	1 x 525 $\mu\text{L}$
R4	1 x 55 $\mu\text{L}$ *
R5	1 x 600 $\mu\text{L}$ *
R6	1 x 50 $\mu\text{L}$
R7	1 x 500 $\mu\text{L}$ *
R8	2 x 1,5 mL
R9	1 x 1 mL

\*Lyophilized reagents. The volumes described above refer to the final volume after the reagents resuspension, as described in the item PROCESS DESCRIPTION, subitem B (Preparation of the Reagents).

**REAGENTS**

**R1. PCR Solution:** Primer, probe, TRIS-HCl.

**R2. Mix Taq:** DNA Polymerase, dNTPs, MgCl<sub>2</sub>.

**R3. Mix Buffer:** TRIS-HCl.

**R4. IC PCR Solution:** Primer, probe, TRIS-HCl.

**R5. Internal Control:** Plasmid, TRIS-HCl.

**R6. Negative Control:** TRIS-HCl.

**R7. Standard A** (2 x 10<sup>5</sup> copies/μL): Plasmid, TRIS-HCl, EDTA.

**R8. Diluent:** TRIS-HCl, EDTA.

**R9. Water:** DNase/RNase free water.

**EQUIPMENTS AND OPERATIONAL INPUTS****Materials in the kit:**

- Reagents described in the previous item
- Usage Instructions (manual)

**Materials needed but not contained in the kit:**

- 1- Real-time PCR thermocycler.
- 2- Laminar flow hood.
- 3- PCR reaction tubes or PCR reaction plate.
- 4- Powder free latex gloves or similar material.
- 5- Microcentrifuge (3,000 - 12,000 rpm).
- 6- Vortex.
- 7- Micropipettes and sterile filter pipette tips (0.5-10μL, 10-100μL, 100-1000μL).
- 8- Nucleic acid extraction kit.

**TRANSPORTATION AND STORAGE CONDITIONS**

The required storage temperature is -20°C (-10 to -30°C).

**After resuspension of lyophilized reagents, the product is stable for 6 months from the date of resuspension.** Avoid freezing and thawing. The kit may be transported between 2-30°C. Protect from light and avoid moisture.

**SPECIAL CARE**

**1- For professional *in vitro* diagnostic use only.**

**2-** Strictly follow the methodology proposed to obtain accurate results.

**3-** Handle and dispose of all biological samples, reagents and materials as it contain infectious agents. Avoid direct contact with biological samples and reagents. Avoid spills and aerosol. Handle and dispose of wastes in accordance with appropriate security procedures.

**4-** It is require skilled personnel for the molecular biology procedures in order to minimize the risk of erroneous results, degradation of nucleic acids contained in the samples or even sample contamination by amplicons.

**5-** It is necessary to have separate areas for the extraction / preparation of reactions and for the amplification / detection of products. Never introduce an amplification product in the area intended for the extraction or preparation of amplification products.

**6-** All samples and reagents should be handled inside a laminar flow hood. It is necessary to use sterile and DNase/RNase free filter pipettes tips.

**7-** Avoid repeated thawing and freezing of the reagents.

**8-** Store the DNA samples at -20°C in case they will not be used immediately.

**9-** Do not use reagents after expiration date.

**10-** We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

**11-** To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website [www.bioclin.com.br](http://www.bioclin.com.br) or upon request by the SAC (Customer Advisory Service) of Quibasa.

**12-** Do not use the product in case of damaged packaging.

**13-** It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

**SAMPLES**

This kit should be used with DNA samples extracted from whole blood, amniotic fluid, lesions or genital / urinary swab. Other types of samples can be used according to medical recommendations or to laboratories.

Samples must be collected according to the laboratory's recommendations for molecular tests and transported and stored as described below<sup>3</sup>:

Amniotic fluid - between 2 and 8 ° C for 1 day.

Other samples - between 2 and 8 ° C for 3 days.

Use DNA samples suitable for PCR amplification with adequate purity and concentration. Avoid repeated freezing and thawing of the samples.

## PROCESS DESCRIPTION

### A. DNA Extraction

The nucleic acids (DNA) must be extracted in accord to the procedure of the chosen extraction kit.

For control of the extraction process, the **Internal Control (R5)** must be prepared and added to the samples during the extraction as described below:

- 1- Add 4µL of the **Internal Control (R5)** to each tube containing the resuspended samples in extraction/lysis buffer.
- 2- Complete the extraction procedure according to the usage instructions of the extraction kit.

**Note: Never add the Internal Control directly to pure biological sample. This may result in degradation of the sample.**

### B. Preparation of the Reagents

1- Before open, spin (*spin* pulse) each of the following reagents: **PCR Solution (R1)**, **IC PCR Solution (R4)**, **Internal Control (R5)** and **Standard A (R7)**.

2- Resuspend the **Mix Taq (R2)** reagent in 525µL of Reagent **Mix Buffer (R3)**.

3- Resuspend the following reagents: **PCR Solution (R1)**, **IC PCR Solution (R4)** and **Internal Control (R5)** with **Water (R9)** reagent, in according to the following table:

Reagent	Volume of Water (R9)
PCR Solution (R1)	55 µL
IC PCR Solution (R4)	55 µL
*Internal Control (R5)	600 µL

\*R5 and R7 reagents contain DNA template. They need to be handle (resuspended) in a separated area from the area designed to prepare

other reagents in order to prevent contamination of other reagents.

### C. Dilution of Quantitative Standard

- 1- Resuspend **Standard A (R7)** with 500  $\mu\text{L}$  of **Diluent (R8)**.
- 2- Select 3 microtubes (not supplied in the kit) suitable for serial dilution of **Standard A (R7)**.
- 3- Pipette 90 $\mu\text{L}$  of the **Diluent (R8)** in each microtube and name them as B, C and D respectively.
- 4- Pipette 10 $\mu\text{L}$  of **Standard A (R7)** in microtube B and mix.
- 5- Replace the tip and pipette 10 $\mu\text{L}$  of microtube B in microtube C and mix.
- 6- Replace the tip and pipette 10 $\mu\text{L}$  of microtube C in microtube D and mix.
- 7- At the end of dilution procedure, it will be available 4 standards with the following concentrations:

Standard A -  $2 \times 10^5$  copies/ $\mu\text{L}$

Standard B -  $2 \times 10^4$  copies/ $\mu\text{L}$

Standard C -  $2 \times 10^3$  copies/ $\mu\text{L}$

Standard D -  $2 \times 10^2$  copies/ $\mu\text{L}$

### D. PCR Preparation

- 1- Prior beginning the assay, select sufficient number of microtubes/wells to be used on the reaction, in accord to the number of samples, controls and quantitative standards to be used.
- 2- Prepare the final PCR solution in accord with reactions number required.

Reactivos	1 Reaction	25 Reactions	50 Reactions	100 Reactions
Mix Taq (R2)	10 $\mu\text{L}$	250 $\mu\text{L}$	500 $\mu\text{L}$	1mL
PCR Solution (R1)	1 $\mu\text{L}$	25 $\mu\text{L}$	50 $\mu\text{L}$	100 $\mu\text{L}$
CI PCR Solution (R4)	1 $\mu\text{L}$	25 $\mu\text{L}$	50 $\mu\text{L}$	100 $\mu\text{L}$
Water (R9)	3 $\mu\text{L}$	75 $\mu\text{L}$	150 $\mu\text{L}$	300 $\mu\text{L}$

To prepare different reaction number must be multiplied the volume of reagents for a reaction by the number of required reactions.

**3-** Pipette 15µL of the final PCR solution in microtubes or wells previously determined.

**4-** Add 5µL of extracted DNA sample (the suggested concentration is 5 ng/µL) or 5µL of Quantitative Standards or 5µL of **Negative Control (R6)** in previously determined tubes.

**5-** Mix well.

**6-** The total volume of the reaction is 20µL. Each PCR run should include the relevant controls: Negative Control, Internal Control and Quantitative Standards.

**7-** Transport the tubes to the thermocycler and go to the section E (Thermocycler Settings for Real Time PCR).

**E. Thermocycler Settings for Real-Time PCR**

Check the operation manual of the Real-time PCR thermocycler to program the experiment.

**1- Set the type of experiment:**

Quantitative Test with Standard Curve.

**2- Set the fluorescent probes:**

	Detector	Quencher
HSV	FAM	NFQ-MGB
Internal Control	VIC	

**Note:** The Negative Control and Quantitative Standards do not present the Internal Control, because it is used for extraction and amplification control of the samples.

**3- Set the quantitative standards:**

Standard A - 2 x 10<sup>5</sup> copies/µL

Standard B - 2 x 10<sup>4</sup> copies/µL

Standard C - 2 x 10<sup>3</sup> copies/µL

Standard D - 2 x 10<sup>2</sup> copies/µL

**4- Set the cycling conditions:**

	Temperature	Time	Cycles
1	95°C	2 minutes	1
2	95°C	10 seconds	50
	60°C	60 seconds	

Set “Data Collection” as “stage 2, step 2 (60@0:60)”.

**F. Result Validation****1- Standard curve**

Standard Curve	Permitted Range	Amplification / Detection
Correlation Coefficient ( $R^2$ )	$0,99 \leq R^2 \leq 1,00$	Valid

If the  $R^2$  value does not fall within the limits of range, the result must be considered invalid and the test must be repeated.

**2- Negative Control**

CT Negative Control		Result	Detection
FAM	VIC		
Undetermined	Undetermined	Negative	Valid

**3- Samples**

HSV		Result	Detection
FAM	VIC		
Determined concentration	$25 \leq CT \leq 31$	Positive	Valid
	$CT < 25$ ou $CT > 31$	Positive	Invalid*
Undetermined concentration	$25 \leq CT \leq 31$	Negative	Valid
	$CT < 25$ ou $CT > 31$	Negative	Invalid*

\*The CT values of Internal Control vary with the process conditions, as the extraction efficiency of the DNA/RNA, the concentration of the samples and the thermocycler settings. Therefore, these conditions must be evaluated when the CT values are not appropriated and, in case it is relevant, the results can be validated.

**Example:** In some cases, the samples with high numbers of copies of DNA / RNA can inhibit the amplification of the Internal Control, resulting in CT value outside the optimal range, this does not invalidate the test results.

If the above requirements are not achieved, the assay must be considered invalid and must be repeated.

**G. Result Interpretation**

This kit is able to detect from 10 to 1.000.000 copies of DNA per reaction. The thermocycler software calculates automatically the samples concentration.

**Example:** If the program shows a concentration of 2.00E+005, this means that the sample concentration is  $2.0 \times 10^5$  copies/ $\mu$ L.

Sample Result in Copies/ $\mu$ L (FAM)	Copies per Reaction
$\geq 1 \times 10^6$	> 1.000.000
$2 \leq \text{Quantity} \leq 9 \times 10^5$	Quantity obtained
< 2	< 10

The non-detection of the pathogen DNA does not exclude the presence of infection when the pathogen titer is below the detection limit of this kit. The results provided by this kit must be interpreted by the medical professional responsible, not being the only parameter for the determination of diagnosis and/or treatment of the patient.

The results obtained must be evaluated considering the clinical data and laboratory tests of the patient.

### **PROCESS LIMITATIONS**

Cross contamination that eventually occurs during the sample collection, processing, transportation and storage may give false results.

### **PRODUCT PERFORMANCE**

#### **QUALITY CONTROL**

#### **Accuracy**

#### **COMPARISON OF METHODS AND METHODOLOGY SPECIFICITY**

The **Bio Gene HSV PCR** test was compared with another test for the quantification of the HSV DNA. 25 analyzes were performed and the results were evaluated. Comparing the two tests, the difference between the results obtained was  $\leq 0.5$  log. This shows that the test has good methodological specificity.

#### **Precision**

#### **REPEATABILITY**

The repeatability was calculated from 10 successive determinations, using 2 positive samples with different concentrations, obtaining the following results:

	<b>Sample 1</b>	<b>Sample 2</b>
<b>Average Concentration (Log)</b>	5,138	3,869
<b>Standard Deviation (Log)</b>	0,054	0,048
<b>Coefficient of Variation (%)</b>	1,068	1,248

**REPRODUCIBILITY**

The reproducibility was calculated from 10 successive determinations for 3 consecutive days, using 1 sample with different concentrations, obtaining the following results:

	Day 1	Day 2	Day 3
<b>Average Concentration (Log)</b>	3,869	3,912	3,878
<b>Standard Deviation (Log)</b>	0,048	0,034	0,043
<b>Coefficient of variation (%)</b>	1,248	0,876	1,108

**Clinical Sensitivity**

The **Bio Gene HSV PCR** showed clinical sensitivity of 99.9% and a clinical specificity of 99.9%.

Method		Clinical Data		Total
Bio Gene HSV PCR	Result	Positive	Negative	
	Positive	25	0	25
	Negative	0	25	25
<b>TOTAL</b>		25	25	50

**Analytical Sensitivity**

It was found in the performance study that the analytical sensitivity of the kit is 0,64 log. This means the assay is capable to detect approximately 4 target molecules in 5µL DNA extraction product, used on the amplification reaction.

**DIAGNOSTIC SIGNIFICANCE**

The Herpes Simplex Virus (HSV), which belongs to the *Herpesviridae* family, causes lesions in the skin and mucous mainly of the face, but can infect other parts of the body such as the genital region.

The HSV-1 often cause lesions on the lips, the inside of the mouth and eyes besides, neurological complications such as encephalitis, have also

been related to HSV-1, transmission occurs by direct contact with infected secretions. The HSV-2 causes lesions, blisters and genital pruritus, is usually transmitted sexually or perinatally. The cross-infection of HSV 1 and 2 can happen if there is oral-genital contact.

The clinical symptoms of primary infection is varied, but in many patients the infection is asymptomatic. The HSV can remain in latent form in the host after primary infection.

Detection of human herpesvirus (1 and 2) by PCR according to the literature, has been considered one of the most sensitive techniques, especially in neonatal disease. Therefore, the quantitative detection of HSV by real time PCR is very important due to its accuracy and agility.

### **BIBLIOGRAPHIC REFERENCES**

1. Aldea, C; Alvarez, C.P; Folgueira, L, et al. Detection, quantification and genotyping of Herpes Simplex Virus in cervicovaginal secretions by real-time PCR: a cross sectional survey. J Clin Microbiol. 2002 Mar;40(3):1060-2.
2. Chudzik, E; Karabin, K; Dzieciotkowski, T, et al. [Development of a novel real-time PCR method for detection of HSV types 1 and 2 DNA using HybProbe chemistry]. Med Dosw Mikrobiol. 2010;62(3):255-62. Polish.
3. CLINICAL AND LABORATORY STANDARDS INSTITUTE (CLSI). Collection, transport, preparation and storage of specimens for molecular methods; approved guideline. CLSI document MM13-A. Pennsylvania, USA: Clinical and Laboratory Standards Institute, 2005.
4. Kessler, H.H; Mühlbauer, G; Rinner, B, et al. Comparison of nested PCR and real time PCR of Herpesvirus infections of central nervous system in HIV patients. J Clin Microbiol. 2000 Jul;38(7):2638-42.

### **CUSTOMER SERVICE**




Customer Advisory Service

Phone: 0800 031 5454

E-mail: [sac@bioclin.com.br](mailto:sac@bioclin.com.br)

ANVISA registration for Bio Gene HSV PCR: 10269360243

## UNIVERSAL SYMBOLOLOGY

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED

**BIO** **GENE**

**Bioclin** · **QUIBASA**

 Rua Teles de Menezes, 92 . Belo Horizonte . MG . Brasil . CEP: 31565-130  
Tel +55 31 3439 5454 . Fax +55 31 3439 5455 . [www.bioclin.com.br](http://www.bioclin.com.br)  
FARM. RESP. Sílvia Wandalsen Arndt - CRF MG 7422  
C.N.P.J.: 19.400.787/0001-07 - Indústria Brasileira