CITRATO

REF A001

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FINALIDADE

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USAGE INSTRUCTIONS

FUNCTION

Anticoagulant solution used in coagulation tests. For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Citrate is used for coagulation studies (prothrombin time, platelet aggregation, fibrinogen determination and other clotting factors) to preserve factors V and VIII. It is the most suitable anticoagulant for tests and therapeutic monitoring of heparin and more appropriate for studies of platelet aggregation.

REAGENTS

Anticoagulant Citrate - Store between 15 and 30°C. Contains: Sodium Citrate 730 mmol/L.

PRESENTATION

Presentation	Volume
1	20 mL
2	250 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Test tubes and other materials needed to collect blood. They can be found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage and transportation temperatures should be between 15 to 30°C. Protect from light and avoid moisture.

SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.
- 2- Water used in material cleaning must be recent and free of contaminants.
- 3- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.
- 4- Separate plasma by centrifuging and observe if it is clear and free of fibrins.
- 5- Do not use phosphate-based detergents to clean glassware.
- 6- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 7- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
- 8- Do not use the product in case of damaged packaging.

9- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

TECHNIQUE

Use 1 drop of anticoagulant for each 3 mL of blood. Smoothly homogenize through inversion. Separate plasma by centrifuging (3.000 rpm), for 5 minutes. Transfer supernatant to a test tube and maintain plasma under refrigeration up to the moment of usage.

PROCEDURE LIMITATIONS

See if the plasma is clear and no traces of fibrin, because the presence of fibrin may be a consequence of deficient homogenization or excess blood related to the volume of anticoagulant used.

NUMBER OF TESTS

A001-1 400 Tests / 50 µL of Reagent
A001-2 5000 Tests / 50 µL of Reagent

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

OBELIS S.A.

Bd. Général Wahis, 53
1030 Brussels, Belgium

CUSTOMER SERVICE

Customer Advisory Service
Phone: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for Citrate kit: 10269360081

Review: August/2014

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED