## HCG STRIP

REF K039

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#### FINALIDADE

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### USAGE INSTRUCTIONS

#### FUNCTION

Immunochromatographic method for quick and qualitative determination of the Human Chorionic Gonadotropin (HCG), through the  $\beta$ hCG fraction. For *in vitro* diagnostic only.

#### PRINCIPLE OF ACTION

##### Methodology: Immunochromatography

The method uses anti- $\beta$ hCG monoclonal antibody that reacts with sample of serum and urine. The mixture moves by a chromatographic membrane by capillary action. The samples with concentration superior to 25 mUI/L of  $\beta$ hCG will form a red line in the area where the anti- $\beta$ hCG is immobilized, indicating that the sample is positive. The reagent mixture will continually be absorbed by the membrane until the area of the control antibody, with the formation of another red line, confirming the test was correctly performed.

#### REAGENTS

**Number 1 - Reactive Strip:** Store between 15 and 30°C.

**Do not freeze.** Contains: Anti- $\beta$ hCG monoclonal antibody conjugate with colloidal gold and Control antibody immobilized on the membrane.

#### PRESENTATION

Packings with 25 / 50 / 100 Reactive Strips.

#### EQUIPMENTS AND OPERATIONAL INPUTS

Watch and stopwatch. They can be found at markets specialized in Clinical Analysis Laboratories.

#### TRANSPORTATION AND STORAGE CONDITIONS

The storage and transport temperature should be between 15 and 30°C. Avoid moisture.

#### SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.
- 2- Strictly follow the methodology proposed to obtain exact results.
- 3- The reactive strips do not present contaminant substances, but once the samples are biological, the safety routine must be followed strictly and the use of disposable gloves is indispensable.
- 4- Use clean and dry recipients to settle the samples.
- 5- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

6- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website [www.bioclin.com.br](http://www.bioclin.com.br) or upon request by the SAC (Customer Advisory Service) of Quibasa.

7- Do not use the product in case of damaged packaging.

8- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

#### SAMPLE

**Serum:** the serum samples are stable for 2 days when kept between 2 and 8°C, and for 03 months at -20°C. Hemolyzed and lipemic samples shall be avoided, because they can cause a wrongly positive result. In the presence of a lifted particle centrifuge the sample before carry out test.

**Urine:** the urine samples are stable for 2 days when kept between 2 and 8°C, and for 03 months at -20°C. Prefer to collect first urine in the morning because it contains a higher concentration of  $\beta$ hCG. In not clear urine or containing erythrocytes must be centrifuge before carry out test. Do not use urine with signs of bacterial infection or with pH higher than 6,0.

#### PROCESS DESCRIPTION


##### TECHNIQUE

##### QUALITATIVE PROVE


- 1- The sample must be in a temperature between 15 and 30°C before test begins.
- 2- Take off the reactive strip of the protection case and identify it in an adequate form.
- 3- Place the strip absorbent tip in contact with the sample (serum or urine) during 10 seconds. Do not dip the strip in sample beyond the indicative sample level signs. As an alternative to the immersion of the reactive strip in sample it could pipette 100  $\mu$ L of sample (serum or urine) and disposal over the reactive strip absorbent tip.
- 4- Wait for the formation of lines. The results should be observed after 5 minutes. Do not interpret after 10 minutes.

#### RESULTS


**Positive Test:** formation of two red lines indicating the concentration of  $\beta$ hCG equal or superior to 25 mUI/mL.

**Negative Test:** formation of a single red line in the control region, the negative test or  $\beta$ hCG concentration lower than 25 mUI/mL.

**Invalid Test:** complete absence of line formation or formation of the line only, indicates error in the procedure or deterioration of the reactive strip. In this case, repeat the test using a new reactive strip.

#### PROCEDURE LIMITATIONS

- 1- When used to verify pregnancy, the final diagnosis can not be based on the laboratorial result only, it must correlate the test result with the clinical signs and symptoms.
- 2- The test could be carry out in a probable menstruation day, however very low concentrations of  $\beta$ hCG can generate negative results. In this case a new test should be carried out after two days to confirm the result.
- 3- The  $\beta$ hCG values lower than 50 mUI/mL can generate wrongly positive results (red shade on reactive strip), should be interpreter along with the signs and symptoms of patient and other confirmatory laboratorial tests.
- 4- Low concentrations of  $\beta$ hCG can occur in normal non-pregnant women.
- 5- The patients with non-pregnant trophoblastic diseases, inclusive carcinoma and the hydatidiform mol, can present positive results in the absence of pregnancy.
- 6- Fake negative results can be obtained in ectopic pregnancy cases, eminent abortion, pregnancy toxemia, due to the possible lowering of the  $\beta$ hCG excretion.
- 7- The presence of heterophile antibodies can produce fake positive results. In the presence of positive results without clinical data confirming it, carry out quantitative test.
- 8- The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

**INTERFERENT**

No interference was observed for Acetaminophen up to 20 mg/dL, Acetylsalicylic Acid up to mg/dL, Ascorbic Acid up to 20 mg/dL, Atropine 20 mg/dL, Caffeine 20 mg/dL, Tetracycline up to 20 mg/dL, Ampicillin up to 20 mg/dL, Albumin up to 20 mg/dL, Bilirubin (serum) up to 40 mg/dL, Bilirubin (urine) up to 2 mg/dL, Glucose up to 2 g/dL and Hemoglobin 20 mg/dL.

**INTERNAL QUALITY CONTROL**

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

**PRODUCT PERFORMANCE****QUALITY CONTROL****Accuracy****METHODS COMPARISON, METHODOLOGICAL SPECIFICITY AND SENSITIVITY**

The HCG Strip kit was compared to other methods for  $\beta$ hCG determination. According with the results of 20 clinical samples, the method showed an excellent correlation. With these results we can conclude that the kit shows good methodological specificity.

**Precision****REPEATABILITY**

The repeatability was calculated from 20 successive determinations, using 10 different samples, obtaining the following results:

Samples	N°. of Repetition	Expected Result	Found result
Serum 01	20	100% Positive	100% Positive
Serum 02	20	100% Positive	100% Positive
Serum 03	20	100% Negative	100% Negative
Serum 04	20	100% Positive	100% Positive
Serum 05	20	100% Negative	100% Negative
Urine 01	20	100% Positive	100% Positive
Urine 02	20	100% Negative	100% Negative
Urine 03	20	100% Positive	100% Positive
Urine 04	20	100% Negative	100% Negative
Urine 05	20	100% Positive	100% Positive

**REPRODUCIBILITY**

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 5 different samples, obtaining the following results:

Samples	N°. of days	N°. of Repetition	Expected Result	Found result
Serum 01	3	20	100% Positive	100% Positive
Serum 02	3	20	100% Positive	100% Positive
Serum 03	3	20	100% Negative	100% Negative
Urine 01	3	20	100% Positive	100% Positive
Urine 02	3	20	100% Negative	100% Negative

**Sensitivity**

The test is capable of detecting  $\beta$ hCG concentrations equal or superior to 25 mUI/mL.

The clinical or diagnostic sensitivity and the methodological sensitivity was determined by the sampling of 20 human serum truly positive. The tests showed 100% positive results, no fake negative was found.

**Specificity**

The test uses a specific monoclonal antibody for the sub-unit of  $\beta$  from the hCG. Thus, no cross reaction are not expected with the FSH, LH and TSH hormones in regular concentrations.

**High Dose Hook Effect**

Hook effect wasn't observed with a high concentration of HCG up to 100 UI/mL.

**NUMBER OF TESTS**

25/50/100 tests

**BIBLIOGRAPHIC REFERENCES**

- 1 - BATZER, F. R., Fertility and Sterility, 4:1, 1980.
- 2 - BRAUNSTEIN, G. D.; RASOR, J.; ADLER, D.; DANZER, H.; WADE, M.E., AM. J. Obstet. Gynecol., 126:681, 1976.
- 3 - BRAUNSTEIN, G. D., et al. Ann. Intern. Med., 78:39-45, 1973.
- 4 - LENTON, E. A.; NEAT, L. M.; SULAIMAN, R., Fertility and Sterility, 37:773-778, 1982.
- 5 - BURTIS, Carl A. Tietz Textbook of Clinical Chemistry, 2ª edic., 2114-2119, 1994.

**QUALITY ASSURANCE**

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

**QUIBASA QUÍMICA BÁSICA Ltda**

Rua Teles de Menezes, 92 - Santa Branca  
CEP 31565-130 - Belo Horizonte - MG - Brasil  
Phone.: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455  
E-mail: bioclin@bioclin.com.br  
CNPJ: 19.400.787/0001-07 - Made in Brazil

EC REP **OBELIS S.A.**

Bd. Général Wahis, 53  
1030 Brussels, Belgium

**CUSTOMER SERVICE**

Customer Advisory Service  
Phone.: 0800 0315454  
E-mail: sac@bioclin.com.br

ANVISA registration for HCG Strip kit: 10269360105

Review: April/2018

**UNIVERSAL SYMBOLOGY**

CATALOG NUMBER

MANUFACTURED BY

BATCH CODE

CONTROL

DATE OF MANUFACTURE

POSITIVE CONTROL

USED BY  
(last day of month)

NEGATIVE CONTROL

TEMPERATURE LIMITATION  
(store at)

BIOLOGICAL RISK

CONTAINS SUFFICIENT  
FOR <N> TESTS

INFLAMMABLE

CONSULT INSTRUCTIONS  
FOR USE

CORROSIVE

IN VITRO DIAGNOSTIC DEVICE

POISON

EUROPEAN AUTHORIZED  
REPRESENTATIVE

CE MARK

KEEP AWAY  
FROM SUNLIGHT

DO NOT USE IF  
PACKAGE IS  
DAMAGED